

## SSR-ER Options:

### 12V Portable Power System Option:

To make the D&S Solar Spectrum Reflectometer Model SSR-ER portable, an optional battery pack can be attached. The battery pack comes with its own rechargeable battery and charging system.

The SSR-ER will operate for up to 5 hours continuous when connected to the 12V Battery Pack.



To make the D&S Solar Spectrum Reflectometer Model SSR-ER portable, an optional battery pack can be attached. The battery pack comes with its own rechargeable battery, charging system, and carrying case.

The SSR-ER will operate for up to 5 continuous hours when connected to the 12V Battery Pack.

**Handle option:** For ease of use when the measurement head must be held in place against the sample.

